



DL040 Panther X Follicle Activation Cap High Density 72 LED, 216 Chips, Red +2 Infrared , Scalp Care Serum

Operating manual



Introduction

Prospera Panther X Follicle Activation 72 LED, 216 Chips, Red & Infrared Cap delivers advanced LED light therapy to support a healthier scalp environment and promote fuller, healthier-looking hair. Featuring 72 LED bulbs with 216 LED chips, red and near-infrared light therapy, 3 brightness levels, and 3 selectable light modes, it offers a personalized treatment experience. With a built-in 2000mAh rechargeable battery, it provides convenient, cordless scalp care at home. The included Scalp Care Serum helps nourish the scalp and complements the LED therapy for an enhanced hair care experience.

Product Specifications

- Applicable for supporting a healthier scalp
- Convenient cap design
- Simple Button Controls

•3 Light Modes:

- Steady Red Light
- Pulsing Red Light
- Breathing Red Light

• 3 Brightness Levels: Low, medium, and high

- Battery capacity: 2,000mAh
- Charging time: approximately 3 hours
- Operation time: 1.5-4 hours
- Charging Indicator: Red flashing while charging; green when fully charged

• Wavelengths:

- Red Light: 660nm
- Infrared Light: 850nm
- LEDs: 72*3pcs (red light: 72 chips, near-infrared light: 72*2chips)
- Auto shut-off: 30 minutes of use
- Voltage: 5V/2A
- Power: 10W
- Charging Port: USB-C
- Main unit: 12.2" * 8.3" * 5.1"
- Item weight: 1.1 lbs
- Box: 12.2" * 7.5" * 4.9"
- UPC 850073782217

Content

CONTENT

Safety Warnings	4
Product Illustration	4
Operating the Product.....	5
Consumer Warranty	6
Warranty Card	6

Benefits

Red Light Therapy

- **Helps Support Hair Follicles :** Helps energize the scalp and support hair follicles, creating a healthier environment for hair growth.
- **Promotes Scalp Circulation :** Helps improve scalp microcirculation, allowing the scalp and hair roots to receive better nourishment.
- **Supports Stronger, Fuller-Looking Hair :** With consistent use, red light may help support stronger, healthier-looking hair and improve the appearance of thinning hair.
- **Helps Maintain Scalp Health :** Red light helps soothe the scalp and reduce the appearance of scalp stress, dryness, or discomfort.

Near-Infrared Light Therapy


- **Deeper Scalp Support:** Near-infrared light penetrates deeper into the scalp to provide enhanced support for hair follicles and scalp tissue.
- **Supports Follicle Vitality:** Helps improve the scalp environment and support follicle vitality for healthier-looking hair.
- **Helps Reduce Scalp Stress:** Near-infrared light helps provide a gentle, soothing scalp care experience and supports overall scalp comfort.
- **Works with Red Light for Better Care:** Near-infrared light works together with red light to deliver deeper, more complete LED scalp therapy.

Safety Warnings

Consult your doctor before using this product if you have any lesions, skin conditions, heart issues, pacemakers, or other health-related concerns.

1. This product is intended for adult use only.
2. Do not let babies and children play with the product.
3. Do not use the product on an infant.
4. If you experience discomfort, stop using the product immediately.
5. Do not wash the unit with water or organic solvents, such as paint thinner, as they may damage the product.

Product Illustration



Operating the Product



Step 1: Preparation (Recommended)

Before each treatment, wash your hair and scalp to remove excess oil, dirt, and styling products. Make sure your hair and scalp are completely dry before use.




Step 2: Apply the Scalp Care Serum

Apply an appropriate amount of Scalp Care Serum to the areas of the scalp that require treatment. Use the included dropper to apply the serum evenly if needed.



Step 3: Turn on The Power & Wear the LED Therapy Cap

Press the Light Mode button and Brightness Level button to select your preferred settings. Put on the cap and adjust the rear strap to achieve a comfortable fit for your head size.



Step 4: After Care

After the treatment is complete, remove the cap. Press and hold the both buttons for 2 seconds to turn off the device.



Step 5: Care & Maintenance

This cap is not waterproof and must not be washed with water. If cleaning is required, wipe the surface gently with a soft, dry, or slightly damp cloth.

Recommend usage: daily

To order a new health and beauty product, please visit www.prosperacorp.com.

Use coupon code: PROS\$5 for \$5 off any product above \$30.

Consumer Warranty

Prospera offers a 90-day standard warranty on massage products. Labor and parts are included during the warranty period. Please provide the model number and dated proof of purchase (sales receipt when contacting Prospera for warranty service.

Non-warranty service is provided on a "per incident" basis. The consumer will verify that the product has failed and provide information for servicing the unit. All applicable repairs, parts, shipping, handling, local tax and "per incident" fees will be charged for non-warranty repairs and support calls.

Warranty Card



Warranty Registration Form

Product Name: _____

Product Model: _____ Purchase Price: _____

Date of Purchase: _____ Retailer: _____

First Name: _____ Last Name: _____

Address: _____

City: _____ State: _____ Zip: _____

Phone: _____ Email: _____

Signature and date: _____

Please mail to:
Prospera Corporation
392 West Larch Rd Suite 31
Tracy, CA 95304-1644

Tel: (925)292-5978
Email: sales@prosperacorp.com